

# BlackJack

by HATCH+MED



- ✓ Only magnetic hospital bed cable
- ✓ Universally compatible with all nurse call systems and beds
- ✓ 37 Individually spring loaded, gold plated and magnetically backed pins
- ✓ Dramatically reduce patient falls and service calls immediately

## Magnetic Breakaway Cables

Eliminate repeated damage to 37 pin cables and precious time to replace them. BlackJack is reliable, durable, and reduces bed-related falls. Protect your patients Nursecall connection.



Scan for more  
info & a free  
sample



info@hatchmed.com



844-428-2463

# BlackJack Compatibility

## Wall Side

	<b>CFR2</b> Right angle 2' magnetic cable	<b>CFR6</b> Right angle 6' magnetic cable	<b>XA</b> Magnetic wall brick	<b>XAS</b> Magnetic wall brick (short)
<b>HillRom 37 Pin BIU or ASBC</b>	YES	YES	NO	YES (only vertical)
<b>Vertical 37 Pin Receptacle</b>	YES	YES	YES	YES
<b>Horizontal 37 Pin Receptacle</b>	YES	YES	NO	YES

## Bed Side

	<b>XM2</b> 2' male magnetic cable	<b>XM6</b> 6' male magnetic cable	<b>XF2</b> 2' female magnetic	<b>XF6</b> 6' female magnetic cable
<b>Stryker &amp; Umano</b>	YES	YES	NO	YES (only vertical)
<b>HillRom &amp; Linet</b>	YES	YES	YES	YES

	<b>XMR6</b> Right angle 6' Male magnetic cable	<b>XFR6</b> Right angle 6' Female magnetic cable	<b>XAF</b> Magnetic Cable Adapter
<b>Stryker &amp; Umano</b>	YES	YES	NO
<b>HillRom &amp; Linet</b>	YES	YES	YES



# BlackJack Compatibility

## Cable Pairing

	<b>XM2</b> 2' male magnetic cable	<b>XM6</b> 6' male magnetic cable	<b>XF2</b> 2' female magnetic	<b>XF6</b> 6' female magnetic cable
<b>CFR2</b> Right angle 2' magnetic cable	NO	YES	NO	YES
<b>CFR6</b> Right angle 6' magnetic cable	YES	NO	YES	NO
<b>XA</b> Magnetic wall brick	YES	YES	YES	YES
<b>XAS</b> Magnetic wall brick (short)	YES	YES	YES	YES

	<b>XMR6</b> Right angle 6' Male magnetic cable	<b>XFR6</b> Right angle 6' Female magnetic cable	<b>XAF</b> Magnetic Cable Adapter
<b>CFR2</b> Right angle 2' magnetic cable	NO	NO	YES
<b>CFR6</b> Right angle 6' magnetic cable	NO	NO	YES
<b>XA</b> Magnetic wall brick	YES	YES	YES
<b>XAS</b> Magnetic wall brick (short)	YES	YES	YES



# BlackJack Specifications

## General Info

<b>Cable</b>	<ul style="list-style-type: none"> <li>UL Listed AVLV2 wires</li> <li>37 conductors 24 AWG</li> <li>6.5mm OD PVC jacket</li> </ul>
<b>Connectors</b>	<ul style="list-style-type: none"> <li>DP37 Male &amp; Female</li> <li>Gold plated contacts</li> <li>Spring loaded male connector pins</li> </ul>
<b>Electrical Ratings</b>	<ul style="list-style-type: none"> <li>100 volt-amps, 5 amps at max rated voltage and frequency</li> <li>30VAC 60Hz or 42.4VDC</li> </ul>

## Regulatory

<b>Certifications</b>	<ul style="list-style-type: none"> <li>FDA Registered</li> <li>UL 1069</li> <li>CSA 22.2 No. 205-12</li> </ul>
-----------------------	--



Scan for more info & help picking the best cable for your hospital

