

SnapJack

by HATCH+MED



- ✓ Magnetic breakaway pillow speaker adapter cables are guaranteed to reduce damage and costly repairs to patient station
- ✓ Improves overall patient to nurse communication by reducing risks of unnoticed damage
- ✓ Retrofit and magnetize your existing configuration
- ✓ Compatible with all nurse call systems

Magnetic Breakaway Cables

Eliminate damage to pillow speakers and patient stations with SnapJack. SnapJack provides a reliable, damage-free magnetic connection that protects both the equipment from damage and the patient from broken lines of communication with their nursing team.



Scan for more info & a free sample



info@hatchmed.com

HATCH+MED



844-428-2463